

FMT

Medical Sensors, Cables and Accessories



EN ISO
9001:2008 / 13485:2003

Humidifier

HEATER WIRE ADAPTERS

REV.02 22.05.2015
YMC EDITION



HUMIDIFIER HEATER WIRE ADAPTERS

- HWA Series Heater Wire Adapters are used to provide electrical energy from the humidifier to the heater wire of disposable or reusable heated breathing circuits. Dual or Single Limb Adapters are available. HWA Series Heater Wire Adapters do not have surge and over-current protection circuit.

▼ Breathing Circuit Connector

▼ Device Side Monitor



- REF: **HWA-001** (OEM P/N **F&P 900MR801**, DRÄGER **8414992**)
Heater Wire Adapter for RT Series disposable heated breathing circuits. For use with MR850, HC550 Humidifiers.



- REF: **HWA-002** (OEM P/N **F&P 900MR800**, DRÄGER **MP01006**)
Heater Wire Adapter for RT Series disposable dual heated breathing circuits. For use with MR850, HC550 Humidifiers.



- REF: **HWA-003** (OEM P/N **F&P 900MR858**, DRÄGER **8414968**)
Heater Wire Adapter for reusable heated breathing circuits with LEMO connector. For use with MR850, HC550 Humidifiers.



- REF: **HWA-004** (OEM P/N **F&P 900MR856**)
Heater Wire Adapter for reusable dual heated breathing circuits with LEMO connector. For use with MR850, HC550 Humidifiers.



- REF: **HWA-005** (OEM P/N **F&P 900MR859**, DRÄGER **MP01001**)
Heater Wire Adapter for disposable heated breathing circuits with Molex connector. For use with MR850, HC550 Humidifiers.



- REF: **HWA-006** (OEM P/N **F&P 900MR855**)
Heater Wire Adapter for disposable dual heated breathing circuits Molex connector. For use with MR850, HC550 Humidifiers.

HUMIDIFIER HEATER WIRE ADAPTERS

▼ Breathing Circuit Connector

▼ Device Side Monitor



- REF: **HWA-007** (OEM P/N **F&P 900MR901**)
Heater Wire Adapter for RT Series disposable heated breathing circuits. For use with MR730 Humidifiers.



- REF: **HWA-008** (OEM P/N **F&P 900MR900**)
Heater Wire Adapter for RT Series disposable dual heated breathing circuits. For use with MR730 Humidifiers.



- REF: **HWA-009** (OEM P/N **F&P 900MR558**)
Heater Wire Adapter for reusable heated breathing circuits with LEMO connector. For use with MR730 Humidifiers.



- REF: **HWA-010** (OEM P/N **F&P 900MR556**)
Heater Wire Adapter for reusable dual heated breathing circuits with LEMO connector. For use with MR730 Humidifiers.



- REF: **HWA-011** (OEM P/N **F&P 900MR554**)
Heater Wire Adapter for disposable heated breathing circuits with Molex connector. For use with MR730 Humidifiers.



- REF: **HWA-012** (OEM P/N **F&P 900MR555**)
Heater Wire Adapter for disposable dual heated breathing circuits Molex connector. For use with MR730 Humidifiers.

Heater Wire Breathing Circuit Connectors



RT / F&P Type
(Disposable)




RT / F&P Type
(Disposable)



Molex Type
(Disposable)



LEMO Type
(Reusable)



Exclusively Distributed By:



Saudi Arabia Contact:

ENG. BENJAMIN TED A. SY
Mobile No.: (+966) 55 904 4546
Email Add.: ymc@yalate.com



Metko Medikal ve Tıbbi Cihazlar Dış Ticaret Ltd. Şti.

Headquarter & Factory

İvedik O.S.B. Ağaç İşleri Sanayi Sit. 1354. Cad. 1358.Sok. No:9 06378 Yenimahalle **Ankara / TÜRKİYE**
Tel (pbx) : +90 312 387 12 46 Fax : +90 312 387 12 51

Branch

Atatürk Mah. Sedef Cad. 38 Ada Ata Plaza 3-5 Blokları Kat:2 No:57 34758 Ataşehir **İstanbul / TÜRKİYE**
Tel (pbx) : +90 216 455 91 21 Fax : +90 216 455 48 24