

FMT

Medical Sensors, Cables and Accessories



NIBP CUFFS & HOSES

REV.25 09.03.2016
YMC EDITION



NIBP CUFFS AND HOSES

Reusable

FMT offers a complete range of sizes and configurations in order to meet your specific monitoring needs. Single and double tube cuffs in Thigh, Large Adult, Adult, Small Adult, Pediatric, Small Pediatric and Infant sizes are available. FMT reusable cuffs have both bladderless and with bladder options.

FMT reusable cuffs have clear index lines, range markings and symbols to ensure easy and correct application for maximum measurement accuracy. FMT cuffs immediately and clearly indicates whether or not the correct size being used.




FMT cuffs with bladder are made of tightly woven soft Nylon. Critical points of the cuffs are double stitched for extra durability. Nylon fabric of FMT reusable cuffs with bladder provides easy cleaning and disinfection. Bladders (inflation bags) are available in Non-Latex (PVC and Neoprene).

FMT cuffs with bladder are antimicrobial coated for reducing the risk of cross contamination.

FMT reusable bladderless cuffs are made of soft Nylon which provides easy cleaning, disinfection and durability. FMT reusable bladderless cuffs are Non-Latex.

Velcro® hook and loop fasteners which are precisely positioned and evenly match in holding strength to provide the optimal number of open – close cycles.



All FMT Products meet the requirements of the MDD 93/42/EEC&MDD 2007/47/EC and CE marked.

NIBP CUFFS AND HOSES

Reusable

NIBP Reference Table

Size		Reference No. for Cuffs with Bladder		Reference No. for Bladderless Cuffs		Arm Circumference (cm)
		Single Lumen (Tube)	Double Lumen (Tube)	Single Lumen (Tube)	Double Lumen (Tube)	
Thigh	6	FMT – S6	FMT – D6	FMT – SB6	FMT – DB6	40.6 - 66
Large Adult	5	FMT – S5	FMT – D5	FMT – SB5	FMT – DB5	33 - 47
Adult Extra Long	4L	FMT – S4L	FMT – D4L	FMT – SB4L	FMT – DB4L	25 - 35
Adult	4	FMT – S4	FMT – D4	FMT – SB4	FMT – DB4	25 - 35
Small Adult	3	FMT – S3	FMT – D3	FMT – SB3	FMT – DB3	20 - 26
Pediatric	2	FMT – S2	FMT – D2	FMT – SB2	FMT – DB2	14 - 20
Infant	1	FMT – S1	FMT – D1	FMT – SB1	FMT – DB1	10 - 14

Length of tubes are 300 mm, outer diameter is Ø 7.8 mm and inner diameter is Ø 4.0 mm



FMT Cuff Bladderless Type

FMT Cuff with Bladder Type

Above part numbers are for cuffs without connectors. For the determining of the P/N (or order code) of the cuffs with connector please chose the size, type and tube condition (single or double tube cuff) from the above table then check the manufacturer and the model from the following table. After determining the connector numbers please add connector number(s) at the end of cuff part number.

(Example: FMT-D4/BP1718, FMT-S4/BP12, FMT-DB4/BP1718, FMT-SB4/BP12)

FMT NIBP Hoses

NIBP hoses made of TPU, single and double tube hoses are available. The wide range of available connectors make FMT hoses compatible with most of automated NIBP monitors in the market.

For the determining the P/N (or order code), please refer to the following table. Check the manufacturer model and determine the hose type (single or double hose) and the connector part numbers. Please add connector part numbers at the end of "S" for single hoses and "D" for double hoses. Monitor side connector number(s) first, cuff side connector part numbers second.

(Example: SBP09/BP09, DBP24/BP1818, DBP22/BP1718)



FMT Coiled Extendable NIBP Hoses are made of TPU in 3 meters lenght.

For the determining of the P/N (or order code) of the Coiled Extendable NIBP Hoses please refer to NIBP Cuffs & Hose Connectors Table, check the manufacturer and the model, after determining connector numbers. Please add connector numbers at the end of "CE", the monitor side connector number first, cuff side second.

(Example: CEBP15/BP12)



NIBP CUFFS AND HOSES

Nibp Hoses

Brand / Model	REF / Patient Type	Hose (Monitor Side)	Hose (Cuff Side)	Cuff Connector(s)
Colin BP 10, BP 306, BP 308, BP 408, BP 8800 / 8800C, BP508, BP-S510. BP88NET / NXT / NEXT / SINEXT / SNXT / SNEXT, BP105N / NP / NXMA / NXNE / NXPMA / NXPNE / SN / SNP / SNXMA / SNXNE / SNXPMA / SNXPNE, BX10MA / NE.	SBP32/BP31 Adult / Pediatric	 BP32	 BP31	 BP32
	SBP32/BP30 Infant / Neonate	 BP32	 BP30	 BP20
Datascope 3000, Accutorr 2 / 3 / 4 / Plus Nihon Kohden AP-851PA Module, Lifescope 9, Criticare 507E/P/R, 507N/P, 507S/SD, 508, 5070/P, 1100, 2200 Datascope EL, XL, Passport D, Welch Allyn Atlas 52N00 / 52S00, Fukuda Denshi DS 2000 / 3000 / 5000 / 7100 / 7200 / 7300 / 7001 Mediana M20 M30 YM1000 YM6000	SBP04/BP07 Adult / Pediatric	 BP04	 BP07	 BP04
	SBP04/BP30 Infant / Neonate	 BP04	 BP30	 BP20
GE Datex Cardiocap 5, Cardiocap II, AS/3, CS/3, S/5, S/5 Light,	DBP22/BP1718 Adult / Pediatric	 BP22	 BP1718	 BP1817
	DBP22/BP26 Infant / Neonate	 BP22	 BP26	 BP20
GE Criticon Dinamap+, Invivo MDE Escort Series, M12, M8, Kontron Diascope G2, Kolormon, Micromon, Supermon, Vitalmon, 7150, 7850, Mennen 1100, 2000, 3000, 9000,	DBP0909/BP0909 Adult / Pediatric	 BP0909	 BP0909	 BP0202
	DBP0909/BP26 Infant / Neonate	 BP0909	 BP26	 BP2020
Dräger/Siemens SC6000 / 7000 / 8000 / 9000, Infinity, Delta / DeltaXL, Gamma XL / Gamma XXL, Kappa / Kappa XLT, M540, Vista / Vista XL, Axion Sensis XP	SBP40/BP15 Adult / Pediatric	 BP40	 BP15	 BP12
	SBP40/BP30 Infant / Neonate	 BP40	 BP30	 BP20

NIBP CUFFS AND HOSES

Nibp Hoses

Brand / Model	REF / Patient Type	Hose (Monitor Side)	Hose (Cuff Side)	Cuff Connector(s)
Drager - Siemens SC6000-7000/8000 9000 Infinity	SBP40/BP46 Adult / Pediatric	 BP40	 BP46	 BP45
	DBP24/BP1717 Adult / Pediatric	 BP24	 BP1717	 BP1818
	Not Available Adult / Pediatric	 BP24	N/A	 BP44
GE Tram 400 / 450 / 451 / 800 / 850 / 851, Coro 120, Marquette B650, B850, Dash 2000 / 2500 / 3000 / 4000 / 5000, Eagle 1000, Solar 8000 TRAM	DBP54/BP26 Infant / Neonate	 BP54	 BP26	 BP2020
	DBP25/BP0909 Adult / Pediatric	 BP25	 BP0909	 BP0202
	DBP25/BP26 Infant / Neonate	 BP25	 BP26	 BP2020
GE Marquette Tram 100 / 200 / 250 / 300 / 350, Solar, Eagle 4000, Corometrics 120	DBP22/BP1717 Adult / Pediatric	 BP22	 BP1717	 BP1818
	DBP22/BP26 Infant / Neonate	 BP22	 BP26	 BP2020
MEK MP700, MP1000 Schiller Argus VCM	SBP09/BP09 Adult / Pediatric	 BP09	 BP09	 BP02
	SBP09/BP30 Infant / Neonate	 BP09	 BP30	 BP20

NIBP CUFFS AND HOSES

Nihp Hoses

Brand / Model	REF / Patient Type	Hose (Monitor Side)	Hose (Cuff Side)	Cuff Connector(s)
Nihon Kohden BSM2300 Lifescope i/L BSM4100 Lifescope P BSM5100 Lifescope A BSM9510 Lifescope M BSM9880 Lifescope S	DBP35/BP18 OEM P/N: YN901P Adult / Pediatric	 BP35	 BP18	 BP17
	DBP36/BP30 OEM P/N: YN921P Infant / Neonatal	 BP36	 BP30	 BP20
Nihon Kohden OPV-1500 Lifescope N	DBP37/BP18 OEM P/N: YN101T Adult / Pediatric	 BP37	 BP18	 BP17
	DBP37/BP30 OEM P/N: YN121T Infant / Neonatal	 BP37	 BP30	 BP20
Philips All except A1/A3, Biolight AnyView2E, AnyView3, AnyView5, AnyView6, AnyView8, EMS, M69, Q3, Q5, Q7, V6, Goldway G60, G70, G80, Mindray BeneViewT1 / T5 / T6 / T8, iMEC 8-10-12, iPM 8-10-12, iPM9800, PM6000, VS600, VS900, BeneHeart D6, MPM Module, DPM 6/7 Üzümcü Visio Mod	SBP16/BP15 Adult / Pediatric	 BP16	 BP15	 BP12
	SBP16/BP30 Infant / Neonatal	 BP16	 BP30	 BP20
Adecon DK 7000, DK 8000/D/L/M/N/S Philips Compact A1/A3, Nellcor 3100 / 4000 Datascope , Spacelabs Ultraview SL2200 / 2400 / 2600 / 2700 / 2800, Elance Elite VSM series 5 / 7, mCAR 300 VSM, Scout, Xprezzon, Charmcare Prizm 3, Kontron , Mindray M7000 / 8000 / 8000A / 8500 / 9000 / 9000A / 9500, DPM 3/4/5 Biosys BPM 200/300/700, Guardian, Sentry-NS Bionet BM3 / 5 / 7 Edan M3/3A/3B, iM50/M50, M8, M8A, M8B, iM80/M80, M9, iM9A/M9A, Goldway G30 / 40, UT4000A / F, UT6000A Mennen Envoy, Menntor X7, Vitalogik 4000 / 4500 / 5000 / 5500 / 6000, MennMove 3/50/70 Criticare 506DN, 506N3 VitalCare Series, 8100E nGenuity, 8100E1 nGenuity, 8100H nCompass, 8500 Poet IQ, CSI 507 Zoll M Series Üzümcü Visio Infunix IP1020 / 3010 / 3050 / 4050 Innomed InnoCare T-12, CardioAid 360-B Spengler Monica 200/2000 Contec CMS 5000 / 5000C / 5100 / 6000 / 6800 / 7000 / 7000Plus / 8000 / 9000, HMS 6500 Cardiac Science 5300 / 5350 Petas KMA 175/275/375/800, PM 100/150 Emtel FX 2000P / 3000 / 3000MD Schiller Argus LCM/LCM Plus	SBP15/BP15 Adult / Pediatric	 BP15	 BP15	 BP12
	SBP15/BP30 Infant / Neonatal	 BP15	 BP30	 BP20
Physio Control Lifepack 20 Cas Medical PPM3, CAS 750, CAS 740, RGB Medical Omicrom FT Surveyor, Omicrom Premium, Omicrom Vision, Dinapro, Oxydina	SBP12/BP15 Adult / Pediatric	 BP12	 BP15	 BP12
	SBP12/BP30 Infant / Neonatal	 BP12	 BP30	 BP20

NIBP CUFFS AND HOSES

Nibp Hoses

Brand / Model	REF / Patient Type	Hose (Monitor Side)	Hose (Cuff Side)	Cuff Connector(s)
Schiller Argus LCM, Argus LCM +, Argus PRO, Argus LCX	SBP15/BP12 Adult / Pediatric	 BP15	 BP12	 BP15
	SBP15/BP30 Infant / Neonate	 BP15	 BP30	 BP20
Welch Allyn Atlas	SBP04/BP09 Adult / Pediatric	 BP04	 BP09	 BP02
	SBP04/BP30 Infant / Neonate	 BP04	 BP30	 BP20
Comen C20 / 30 / 50 / 60 / 70 / 80 / 100, Infinium Omni 2 / 3 / Express / Express VS Petas KMA 900 Heal Force Deluxe 70/100/140, Advance 120/150 UP800, PC900	SBP43/BP15 Adult / Pediatric	 BP43	 BP15	 BP12
	SBP43/BP30 Infant / Neonate	 BP43	 BP30	 BP20
Dixtal DX2010 / DX2020 / DX2021 / DX2022 / DX2023 / DX2025 / DX2405 / DX7100	DBP22/BP4242 Adult / Pediatric Infant / Neonate	 BP22	 BP4242	 BP4141

NIBP CUFFS AND HOSES

Nibp Connectors



• BP 00*



• BP 01**



• BP 02 • BP 09



• BP 03 • BP 04 • BP 06



• BP 04 • BP 03 • BP 07



• BP 06 • BP 07 • BP 03



• BP 07 • BP 04



• BP 09 • BP 09-P • BP 02



• BP 09-P • BP 09 • BP 02



• BP 12 • BP 12-P • BP 15



• BP 12-P • BP 12 • BP 15



• BP 15 • BP 12 • BP 12-P




• BP 16



• BP 17 • BP 18



• BP 18 • BP 17



• BP 20 • BP 23 • BP 30



• BP 21 • BP 21-P • BP 21



• BP 22



• BP 23 • BP 20 • BP 30

* Please note that BP00 connector is only used for closing cuff tubes during cleaning.

** BP01 is used for closing cuff tube with BP12 connector during cleaning.

• Equivalent
• Mate

NIBP CUFFS AND HOSES

Nibp Connectors



• BP 24



• BP 25



• BP 26 • BP 20



• BP 30 • BP 20



• BP 31 • BP 32



• BP 32 • BP 31



• BP 35



• BP 36



• BP 37



• BP 40



• BP 41



• BP 42



• BP 43



• BP 44



• BP 45



• BP 46



• BP 50



• BP 54



• Y-ADP



• BP 18Y



• BP 30Y

• Equivalent
• Mate

NIBP CUFFS AND HOSES

Nibp Y Adapters



Y-ADP BP12/BP1718



Y-ADP BP02/BP1718



Y-ADP BP18/BP1718



Y-ADP BP04/BP1718



Y-ADP BP12/BP1717



Y-ADP BP02/BP1717



Y-ADP BP18/BP1717



Y-ADP BP04/BP1717



Y-ADP BP12/BP0909



Y-ADP BP02/BP0909



Y-ADP BP18/BP0909



Y-ADP BP04/BP0909



Y-ADP BP12/BP3030



Y-ADP BP02/BP3030



Y-ADP BP18/BP3030



Y-ADP BP04/BP3030



Y-ADP BP30/BP0909



Y-ADP BP30/BP1718



Y-ADP BP30/BP0202



Y-ADP BP30/BP1818



Y-ADP BP15/BP0202



Y-ADP BP15/BP1718



Y-ADP BP15/BP1818



Y-ADP BP18/BP1818

NIBP CUFFS AND HOSES

Single Patient Use

■ Features

- Designed for single patient use to eliminate cross-contamination.
- Durable hook-and-loop closure allows repeated repositioning for ease of use.
- Extra soft and comfortable TPU material is durable and 100% latex-free.
- Wide range of sizes to accommodate all patients. Cuffs are available in 12 sizes.
- In Neonate sizes color coded printing on each cuff makes it easy to identify cuff size.
- Dual or single tube styles to suit the configuration of your blood pressure equipment.
- Range guide and index marking on the cuff help to determine the proper size cuff for the patient.
- Adapters available for any monitor on the market.
- Cuffs have 20 cm tube/s.
- Packaged 10 or 40 per case.



■ Size and compatibility

#1 to #5 sizes Neonatal cuffs are supplied with a Luer slip connector. Standard male Luer connectors allow "FMT neonatal" blood pressure cuffs to be adapted to a wide range of most major monitor manufacturers.

#6 to #12 sizes cuffs are shipped without connector or with the plastic connectors attached that allow them to be used with virtually any NIBP monitor or sphygmomanometer.

#6 to #12 sizes Single Patient Use blood pressure cuffs are available with any of the following connectors in either one- or two-tube versions.



"FMT Single Patient Use blood pressure cuffs are designed to reduce the spread of infection through cross-contamination. FMT offers a complete line of single patient use NIBP cuffs from Large Adult size to Neonatal small size. All cuffs are manufactured from a soft, durable, white TPU material which can be quickly and easily wiped clean. The cuff has been manufactured with rounded corners to minimize discomfort."



NIBP CUFFS AND HOSES

Single Patient Use

FMT NeoCuff Single Patient Use Neonatal NIBP Cuffs – Single Tube

REF:	Size	Cuff Dimensions L x W (cm)	Cuff Range (cm)	Color
FMT-DS5	#5 Infant	19.1 x 5.6	8.0 – 15.0	Orange
FMT-DS4	#4 Neonatal Extra Large	18.1 x 4.8	7.0 – 13.0	Red
FMT-DS3	#3 Neonatal Large	14.6 x 4.1	6.0 – 11.0	Green
FMT-DS2	#2 Neonatal Medium	11.9 x 3.1	4.0 – 8.0	Light Blue
FMT-DS1	#1 Neonatal Small	10.0 x 2.5	3.0 – 6.0	Navy Blue

FMT NeoCuff Single Patient Use Neonatal NIBP Cuffs – Double Tube

REF:	Size	Cuff Dimensions L x W (cm)	Cuff Range (cm)	Color
FMT-DD5	#5 Infant	19.1 x 5.6	8.0 – 15.0	Orange
FMT-DD4	#4 Neonatal Extra Large	18.1 x 4.8	7.0 – 13.0	Red
FMT-DD3	#3 Neonatal Large	14.6 x 4.1	6.0 – 11.0	Green
FMT-DD2	#2 Neonatal Medium	11.9 x 3.1	4.0 – 8.0	Light Blue
FMT-DD1	#1 Neonatal Small	10.0 x 2.5	3.0 – 6.0	Navy Blue

FMT Single Patient Use NIBP Cuffs – Single Tube

REF:	Size	Cuff Dimensions L x W (cm)	Cuff Range (cm)	Color
FMT-DS12	#12 Large Adult	63.0 x 17.0	33.0 – 47.0	White
FMT-DS11	#11 Adult	51.0 x 14.0	25.3 – 34.3	White
FMT-DS10	#10 Small Adult	40.0 x 11.4	20.0 – 27.0	White
FMT-DS9	#9 Pediatric	38.0 x 8.6	15.8 – 21.3	White
FMT-DS8	#8 Small Pediatric	30.0 x 7.2	12.4 – 16.8	White
FMT-DS7	#7 Infant	25.0 x 6.0	9.8 – 13.3	White
FMT-DS6	#6 Neonate	20.0 x 5.0	7.7 – 10.5	White

FMT Single Patient Use NIBP Cuffs – Double Tube

REF:	Size	Cuff Dimensions L x W (cm)	Cuff Range (cm)	Color
FMT-DD12	#12 Large Adult	63.0 x 17.0	33.0 – 47.0	White
FMT-DD11	#11 Adult	51.0 x 14.0	25.3 – 34.3	White
FMT-DD10	#10 Small Adult	40.0 x 11.4	20.0 – 27.0	White
FMT-DD9	#9 Pediatric	38.0 x 8.6	15.8 – 21.3	White
FMT-DD8	#8 Small Pediatric	30.0 x 7.2	12.4 – 16.8	White
FMT-DD7	#7 Infant	25.0 x 6.0	9.8 – 13.3	White
FMT-DD6	#6 Neonate	20.0 x 5.0	7.7 – 10.5	White



Single Tube Cuffs



BP12-P

CAS 740-750, Criticare 508, Datascope Accutor Plus & Passport 2, Philips/ HP, Siemens, Spacelabs.



BP03

BCI, Colin, Datascope, Passport XG and Expert, Nihon Kohden, Lifescope 9&12, Welch Allyn Vital Signs, Atlas



BP02

Welch Allyn Vital Signs, Atlas



BP18

Welch Allyn Protocol Propaq



Double Tube Cuffs



BP0202

Critikon, Dinamap, Criticare 507, Invivo Research, MDE, Nihon Kohden Lifescope L, Spacelabs



BP1718


Datex Ohmeda



BP1818

GE Medical Eagle 4000, Tram 100/200/300

Please add the connector numbers at the end of the cuff reference number for the complete cuff part number. Example: FMT-DS12/BP12, FMT-DD10/BP0202.



Exclusively Distributed By:



Saudi Arabia Contact:

ENG. BENJAMIN TED A. SY
Mobile No.: (+966) 55 904 4546
Email Add.: ymc@yalate.com



Metko Medikal ve Tıbbi Cihazlar Dış Ticaret Ltd. Şti.

Headquarter & Factory

İvedik O.S.B. Ağaç İşleri Sanayi Sit. 1354. Cad. 1358.Sok. No:9 06378 Yenimahalle **Ankara / TÜRKİYE**
Tel (pbx) : +90 312 387 12 46 Fax : +90 312 387 12 51

Branch

Atatürk Mah. Sedef Cad. 38 Ada Ata Plaza 3-5 Blokları Kat:2 No:57 34758 Ataşehir **İstanbul / TÜRKİYE**
Tel (pbx) : +90 216 455 91 21 Fax : +90 216 455 48 24